EN: Indication for use

Media for immobilization and isolation of sperm prior to intracytoplasmic sperm injection, ICSI.

Product Description

ICSI™ is a viscous sperm handling solution containing polyvinylpyrrolidone (PVP) and recombinant human albumin (rHA).

ICSI™ is antibiotic free.

For use after equilibration at +20 ±5°C and ambient atmosphere.

Storage instructions and stability

Store dark at +2 to +8°C.

ICSI™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The following is the general procedure for using ICSI™. For information on recommended use of ICSI™, please see "G Series™ Manual" or "Vitrolife Fertility System Manual". The manuals are available on request from Vitrolife, and may be ordered through the company's internet website, www.vitrolife.com. Direction for use in other languages can also be found on the website.

Prepare warmed and equilibrated dishes with droplets of handling-medium and a droplet of ICSI™ in the center.

Cover with OVOIL™.

Place the oocytes into the handling medium droplets, one oocyte per droplet.

Place a small volume (1–2 µL) of prepared sperm suspension into the center of the ICSI™-droplet. Incubate the dishes for a few minutes to allow the sperm to migrate to the outer perimeter of the droplet.

Prime the injection pipette with ICSI™.

Immobilize individual sperm by using the injection pipette to crush the membrane of the sperm tail. Aspirate the immobilized sperm.

Move the oocyte droplet into the field of view. Use a holding pipette to secure the oocyte for injection. Inject the sperm with the least amount of ICSI™.

Place the injected oocytes in culture medium.

Specifications

Sterile filtered	SAL 10 ³
Bacterial endotoxins (LAL assay) [EU/mL]	< 0,25
Human sperm motility recovery assay	≥40% of initial motility
Microscopic evaluation	Passes test
Mouse Embryo Assay	Not performed
LOT specific test results are available on the Certificate of Analysis provided with each delivery.	

Symbols

	Sterilized using aseptic processing techniques
	Temperature limitation
	Do not re-use, discard after procedure.
	Use by - see label.
	Caution: Consult accompanying documents.
	Catalog number
	Batch code

Precautions

Do not use ICSI™ if it appears cloudy.

ICSI™ contains polyvinylpyrrolidone and recombinant human albumin, but no antibiotics.

Vitrolife recommends that media be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

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Vitrolife Sweden AB, Box 9080 (Gustaf Werners gata 2), SE-400 92 Göteborg, Sweden.
Vitrolife Inc., 3601 South Inca Street, Englewood CO 80110, USA.
Vitrolife Pty. Ltd., 55 Flemington Rd, North Melbourne VIC 3051, Australia.
Vitrolife Sweden AB